


A large, close-up photograph of an ENFit enteral syringe. The syringe is filled with a bright orange liquid. The plunger is visible at the top, and the barrel has a black ring near the top. The background is a soft, out-of-focus white.

Baxter

EVERY
CONNECTION
COUNTS



ENFit® Enteral Syringes

SYRINGES FOR A NEW ERA IN ENTERAL ADMINISTRATION

For 45 years, Baxter has offered a comprehensive line of syringes designed to **prevent misconnection** with other devices. In an effort to prevent wrong route delivery of fluids and gases, the International Organization for Standardization (ISO), in collaboration with the FDA, AAMI, and clinicians, has developed a new standard to reduce the likelihood of tubing misconnections, ISO 80369. Baxter, in support of this global initiative to **enhance patient safety**, is introducing a complete line of ENFit enteral syringes and accessories that complies with ISO 80369-3 specific for enteral devices.




COMPLETE LINE OF ENFIT SYRINGES

Sizes available from **60 mL to 0.5 mL** to support enteral-route medication administration. The clear barrel, dark print, and ENFit design allow for easier drawing of medication and accurate dose delivery. **For use with standard liquids or viscous fluids.** Syringes shown actual size.

UNIQUE TIP DESIGN

The new ENFit design allows the clinician to **feel safer** when using the syringe for enteral administration.

- ISO 80369-3 compliant; **designed to prevent misconnections**
- Connector design provides **secure connection** to tube set
- Low Dose Tip (LDT) reduces the dead space, ensuring **delivery accuracy** for small doses (available in 6 mL and smaller)



PHARMACY BOTTLE CAPS

Baxter's line of Pharmacy Bottle Caps facilitate easier and more efficient filling of syringes from a medication bottle. Five of the most popular sizes are designed specifically for the Baxter ENFit syringes.



TIP CAPS

Unique ENFit tip caps designed to fit specifically to the Baxter enteral syringes. Securely cap the ENFit enteral syringe after filling for transportation and medication administration.



CARE FOR THE SMALLEST PATIENTS

ENFit syringes that comply with the ISO standard and are suitable for use with **neonatal, pediatric, and adult** patient populations.

SPECIFICATIONS

ISO 80369-3 COMPLIANT

SYRINGE MATERIALS

- polypropylene
- polydimethylsiloxane
- silicone
- black ink

TIP CAP MATERIALS

- polypropylene
- white colorant

PHARMACY BOTTLE CAP MATERIALS

- polypropylene
- purple colorant

ORDERING INFORMATION



REF

0.5 mL Non-Sterile ENFit Syringe, Orange	500	H93899005
1 mL Non-Sterile ENFit Syringe, Orange	500	H93899101
3 mL Non-Sterile ENFit Syringe, Orange	500	H93899103
6 mL Non-Sterile ENFit Syringe, Orange	500	H93899106
12 mL Non-Sterile ENFit Syringe, Orange	500	H93899112
20 mL Non-Sterile ENFit Syringe, Orange	200	H93899120
60 mL Non-Sterile ENFit Syringe, Orange	200	H93899160
Pharmacy Bottle Cap B, 20 mm, Purple	25	H9386102EN
Pharmacy Bottle Cap C, 22 mm, Purple	25	H9386103EN
Pharmacy Bottle Cap D, 24 mm, Purple	25	H9386104EN
Pharmacy Bottle Cap E, 28 mm, Purple	25	H9386105EN
Pharmacy Bottle Cap I, 38 mm, Purple	25	H9386109EN
Tip Cap, White	600	H93853305

MAKE SAFE ORAL CONNECTIONS

We've created an Oral Liquid Medication Delivery portfolio with a full range of syringes, sizes, packaging configurations and accessories so that you can decide the drug delivery system that works best for your pharmacy. Our goal is to help you streamline your pharmacy practices. That is why we developed a complete system for preparing and administering oral medications with safety and efficiency in mind.



For more information about ISO 80369-3 and the Stay Connected Initiative, visit www.stayconnected.org or contact Baxter at xxx.xxx.xxxx

Baxter is a trademarks of Baxter International Inc., or its subsidiaries.
ENFit is a registered trademark of GEDSA used with their permission.

Baxter